

## Anti-UPP1 Rabbit Monoclonal Antibody

Catalog Number: M10128

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




## Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-UPP1 Rabbit Monoclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                               |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Description          | Boster Bio Anti-UPP1 Rabbit Monoclonal Antibody catalog # M10128. Tested in WB application. This antibody reacts with Human, Mouse, Rat.                                                                                                                                                                                                                                                                                                           |
| Application          | WB                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Clonality            | Monoclonal 26U88                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Formulation          | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                                                         |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Uniprot ID           | Q16831                                                                                                                                                                                                                                                                                                                                                                                                                                             |

## Technical Details

|                     |                                               |
|---------------------|-----------------------------------------------|
| Immunogen           | A synthesized peptide derived from human UPP1 |
| Isotype             | IgG                                           |
| Form                | Liquid                                        |
| Concentration       | 0.5mg/ml                                      |
| Purification        | Affinity-chromatography                       |
| Suggested Dilutions | WB 1:500-2000                                 |

## Anti-UPP1 Rabbit Monoclonal Antibody (M10128) Images



Western blot analysis of UPP1 using anti-UPP1 antibody (M10128). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human Hela whole cell lysates, Lane 2: human 293T whole cell lysates, Lane 3: human Jurkat whole cell lysates, Lane 4: human MCF-7 whole cell lysates, Lane 5: rat heart tissue lysates, Lane 6: rat liver tissue lysates, Lane 7: mouse heart tissue lysates, Lane 8: mouse liver tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-UPP1 antigen affinity purified monoclonal antibody (M10128) at 1:500 overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:500 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for UPP1 at approximately 28 kDa. The expected band size for UPP1 is at 34 kDa.

### Anti-UPP1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.