

## Anti-TRIM38 Mouse Monoclonal Antibody [Clone ID: OTI15A3]

Catalog Number: M10133-1


### Overview

Product Name	Anti-TRIM38 Mouse Monoclonal Antibody [Clone ID: OTI15A3]
Reactive Species	Dog, Human, Monkey, Mouse, Rat
Description	Boster Bio TRIM38 mouse monoclonal antibody, clone OTI15A3 (formerly 15A3). Catalog# M10133-1. Tested in IF, IHC, WB. This antibody reacts with Human, Monkey, Mouse, Rat, Dog.
Application	IF, IHC, WB
Clonality	Monoclonal OTI15A3
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	O00635


### Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 1-265 of human TRIM38 (NP_006346) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:200 IHC 1:150 IF 1:100


## Anti-TRIM38 Mouse Monoclonal Antibody [Clone ID: OTI15A3] (M10133-1) Images




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-TRIM38 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).




Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-TRIM38 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Immunohistochemical staining of paraffin-embedded Carcinoma of Human thyroid tissue using anti-TRIM38 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Immunohistochemical staining of paraffin-embedded Human endometrium tissue within the normal limits using anti-TRIM38 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Anti-TRIM38 mouse monoclonal antibody (M10133-1) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY TRIM38.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRIM38 Mouse Monoclonal Antibody [Clone ID: OTI15A3]

For Research Use Only. Not for use in diagnostic procedures.