

## Anti-PBLD Mouse Monoclonal Antibody [Clone ID: OTI4F6]

Catalog Number: M10149-1


### Overview

Product Name	Anti-PBLD Mouse Monoclonal Antibody [Clone ID: OTI4F6]
Reactive Species	Human
Description	Boster Bio PBLD mouse monoclonal antibody, clone OTI4F6 (formerly 4F6). Catalog# M10149-1. Tested in IHC, WB. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal OTI4F6
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P30039


### Technical Details

Immunogen	Full length human recombinant protein of human PBLD (NP_071412) produced in HEK293T cell.
Isotype	IgG2a
Concentration	0.65 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IHC: 1:150

## Anti-PBLD Mouse Monoclonal Antibody [Clone ID: OTI4F6] (M10149-1) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PBLD (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PBLD.



Immunohistochemical staining of paraffin-embedded Human liver tissue within the normal limits using anti-PBLD mouse monoclonal antibody. (M10149-1)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PBLD Mouse Monoclonal Antibody [Clone ID: OTI4F6]

For Research Use Only. Not for use in diagnostic procedures.