

Anti-SGK196 (POMK) Mouse Monoclonal Antibody [Clone ID: OTI2B10]

Catalog Number: M10251


Overview

| | |
|----------------------|--|
| Product Name | Anti-SGK196 (POMK) Mouse Monoclonal Antibody [Clone ID: OTI2B10] |
| Reactive Species | Human |
| Description | Boster Bio SGK196 mouse monoclonal antibody, clone OTI2B10 (formerly 2B10). Catalog# M10251. Tested in IHC, WB. This antibody reacts with Human. |
| Application | IHC, WB |
| Clonality | Monoclonal OTI2B10 |
| Formulation | PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C as received. |
| Host | Mouse |
| Uniprot ID | Q9H5K3 |


Technical Details

| | |
|---------------------|---|
| Immunogen | Human recombinant protein fragment corresponding to amino acids 141-350 of human SGK196 (NP_115613) produced in E.coli. |
| Isotype | IgG1 |
| Concentration | 1 mg/ml |
| Purification | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Suggested Dilutions | WB 1:2000 IHC 1:150 |

Anti-SGK196 (POMK) Mouse Monoclonal Antibody [Clone ID: OTI2B10] (M10251) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SGK196 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SGK196.



Immunohistochemical staining of paraffin-embedded Carcinoma of Human pancreas tissue using anti-SGK196 Mouse monoclonal antibody. (M10251; heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SGK196 (POMK) Mouse Monoclonal Antibody [Clone ID: OTI2B10]

For Research Use Only. Not for use in diagnostic procedures.