

Anti-DPPA4 Antibody (Center)

Catalog Number: M10283

About DPPA4

May be involved in the maintenance of active epigenetic status of target genes. May inhibit differentiation of embryonic cells into a primitive ectoderm lineage.


Overview

Product Name	Anti-DPPA4 Antibody (Center)
Reactive Species	Human
Description	Boster Bio Anti-DPPA4 Antibody (Center) (Catalog # M10283). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q7L190

Technical Details

Immunogen	This DPPA4 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 107-142 amino acids from the Central region of human DPPA4.
Predicted Reactive Species	Bovine
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

Anti-DPPA4 Antibody (Center) (M10283) Images



Anti-DPPA4 Antibody (Center) at 1:1000 dilution + NCCIT whole cell lysates
Lysates/proteins at 20 μ g per lane.
Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/10000 dilution
Predicted band size : 34 kDa
Blocking/Dilution buffer: 5% NFD/MTBST.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DPPA4 Antibody (Center)

For Research Use Only. Not for use in diagnostic procedures.