

Anti-CTDSP2 Mouse Monoclonal Antibody [Clone ID: OTI7F5]

Catalog Number: M10355


Overview

| | |
|----------------------|--|
| Product Name | Anti-CTDSP2 Mouse Monoclonal Antibody [Clone ID: OTI7F5] |
| Reactive Species | Human, Mouse |
| Description | Boster Bio CTDSP2 mouse monoclonal antibody, clone OTI7F5 (formerly 7F5). Catalog# M10355. Tested in IHC, WB. This antibody reacts with Human, Mouse. |
| Application | IHC, WB |
| Clonality | Monoclonal OTI7F5 |
| Formulation | PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C as received. |
| Host | Mouse |
| Uniprot ID | O14595 |


Technical Details

| | |
|---------------------|---|
| Immunogen | Full length human recombinant protein of human CTDSP2 (NP_005721) produced in E.coli. |
| Isotype | IgG1 |
| Concentration | 1 mg/ml |
| Purification | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Suggested Dilutions | WB 1:2000 IHC 1:150 |

Anti-CTDSP2 Mouse Monoclonal Antibody [Clone ID: OTI7F5] (M10355) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CTDSP2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CTDSP2 (1:2000).



Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-CTDSP2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120Å°C for 3min)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CTDSP2 Mouse Monoclonal Antibody [Clone ID: OTI7F5]

For Research Use Only. Not for use in diagnostic procedures.