

## Anti-C1orf19 (TSEN15) Mouse Monoclonal Antibody [Clone ID: OTI5A3]

Catalog Number: M10736


### Overview

Product Name	Anti-C1orf19 (TSEN15) Mouse Monoclonal Antibody [Clone ID: OTI5A3]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio TSEN15 mouse monoclonal antibody, clone OTI5A3. Catalog# M10736. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI5A3
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q8WW01


### Technical Details

Immunogen	Full length human recombinant protein of human TSEN15 (NP_443197) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:250~500

## Anti-C1orf19 (TSEN15) Mouse Monoclonal Antibody [Clone ID: OTI5A3] (M10736) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TSEN15 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TSEN15 (1:500).



Western blot analysis of extracts (35ug) from 4 different cell lines by using anti-TSEN15 monoclonal antibody (1:250).

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C1orf19 (TSEN15) Mouse Monoclonal Antibody [Clone ID: OTI5A3]

For Research Use Only. Not for use in diagnostic procedures.