

## Anti-P4HA3 Mouse Monoclonal Antibody [Clone ID: OTI9A2]

Catalog Number: M10966


### Overview

Product Name	Anti-P4HA3 Mouse Monoclonal Antibody [Clone ID: OTI9A2]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio P4HA3 mouse monoclonal antibody, clone OTI9A2. Catalog# M10966. Tested in IHC. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Monoclonal OTI9A2
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q7Z4N8

### Technical Details

Immunogen	Full length human recombinant protein of human P4HA3 (NP_878907) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	IHC 1:1000

## Anti-P4HA3 Mouse Monoclonal Antibody [Clone ID: OTI9A2] (M10966) Images



Immunohistochemical staining of paraffin-embedded Human skin tissue within the normal limits using anti-P4HA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-P4HA3 Mouse Monoclonal Antibody [Clone ID: OTI9A2]

For Research Use Only. Not for use in diagnostic procedures.