

Anti-Mitochondrial ribosomal protein L11 (MRPL11) Mouse Monoclonal Antibody [Clone ID: OTI3G2]

Catalog Number: M11059-1


Overview

Product Name	Anti-Mitochondrial ribosomal protein L11 (MRPL11) Mouse Monoclonal Antibody [Clone ID: OTI3G2]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio MRPL11 mouse monoclonal antibody, clone OTI3G2 (formerly 3G2). Catalog# M11059-1. Tested in IF, WB. This antibody reacts with Human, Mouse, Rat.
Application	IF, WB
Clonality	Monoclonal OTI3G2
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9Y3B7


Technical Details

Immunogen	Full length human recombinant protein of human MRPL11 (NP_057134) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:

Anti-Mitochondrial ribosomal protein L11 (MRPL11) Mouse Monoclonal Antibody [Clone ID: OTI3G2] (M11059-1) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MRPL11 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MRPL11.



Anti-MRPL11 mouse monoclonal antibody (M11059-1) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY MRPL11.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mitochondrial ribosomal protein L11 (MRPL11) Mouse Monoclonal Antibody [Clone ID: OTI3G2]