

Anti-GPR83 Mouse Monoclonal Antibody [Clone ID: OTI3E8]

Catalog Number: M11089


Overview

Product Name	Anti-GPR83 Mouse Monoclonal Antibody [Clone ID: OTI3E8]
Reactive Species	Mouse, Rat
Description	Boster Bio GPR83 mouse monoclonal antibody, clone OTI3E8. Catalog# M11089. Tested in WB. This antibody reacts with Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI3E8
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9NYM4

Technical Details

Immunogen	Synthetic peptide
Isotype	IgG2b
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500

Anti-GPR83 Mouse Monoclonal Antibody [Clone ID: OTI3E8] (M11089) Images



Western blot analysis of lysates from 293T cells transfected with (1)pcDNA3.1-mGPR83

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GPR83 Mouse Monoclonal Antibody [Clone ID: OTI3E8]

For Research Use Only. Not for use in diagnostic procedures.