

Anti-RPB11 (POLR2J) Mouse Monoclonal Antibody [Clone ID: OTI1B11]

Catalog Number: M11442


Overview

Product Name	Anti-RPB11 (POLR2J) Mouse Monoclonal Antibody [Clone ID: OTI1B11]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio POLR2J mouse monoclonal antibody, clone OTI1B11. Catalog# M11442. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI1B11
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P52435

Technical Details

Immunogen	Full length human recombinant protein of human POLR2J (NP_006225) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500

Anti-RPB11 (POLR2J) Mouse Monoclonal Antibody [Clone ID: OTI1B11] (M11442) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY POLR2J (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-POLR2J (1:500).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RPB11 (POLR2J) Mouse Monoclonal Antibody [Clone ID: OTI1B11]

For Research Use Only. Not for use in diagnostic procedures.