

Anti-CKMT2 monoclonal Antibody

Catalog Number: M11457


Overview

Product Name	Anti-CKMT2 monoclonal Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-CKMT2 monoclonal Antibody catalog # M11457. Tested in WB applications. This antibody reacts with Human,Rat.
Application	WB
Clonality	Monoclonal 2C1
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P17540

Technical Details

Immunogen	Recombinant full length Human CKMT2.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000

Anti-CKMT2 monoclonal Antibody (M11457) Images



Western blot detection of CKMT2 in Rat Brain lysates using CKMT2 antibody (1:1000 diluted) .

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CKMT2 monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.