

Anti-GBGT1 Antibody (Center)

Catalog Number: M11531

About GBGT1

Catalyzes the formation of some glycolipid via the addition of N-acetylgalactosamine (GalNAc) in alpha-1,3-linkage to some substrate. Glycolipids probably serve for adherence of some pathogens.


Overview

Product Name	Anti-GBGT1 Antibody (Center)
Reactive Species	Human
Description	Boster Bio Anti-GBGT1 Antibody (Center) (Catalog # M11531). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8N5D6

Technical Details

Immunogen	This GBGT1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 154-185 amino acids from the Central region of human GBGT1.
Predicted Reactive Species	Mouse
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:2000

Anti-GBGT1 Antibody (Center) (M11531) Images



All lanes : Anti-GBGT1 Antibody (Center) at 1:2000 dilution
Lane 1: Jurkat whole cell lysate
Lane 2: HeLa whole cell lysate
Lane 3: SK-OV-3 whole cell lysate
Lane 4: A2780 whole cell lysate
Lysates/proteins at 20 µg per lane.
Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 40 kDa
Blocking/Dilution buffer: 5% NFDm/TBST.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GBGT1 Antibody (Center)

For Research Use Only. Not for use in diagnostic procedures.