

Anti-MPPB/PMPCB Antibody Picoband® (monoclonal, 9F13E4) HRP Conjugated

Catalog Number: M11793-HRP

About PMPCB

Mitochondrial-processing peptidase subunit beta is an enzyme that in humans is encoded by the PMPCB gene. This gene is a member of the peptidase M16 family and encodes a protein with a zinc-binding motif. This protein is located in the mitochondrial matrix and catalyzes the cleavage of the leader peptides of precursor proteins newly imported into the mitochondria, though it only functions as part of a heterodimeric complex.

Overview

| | |
|----------------------|---|
| Product Name | Anti-MPPB/PMPCB Antibody Picoband® (monoclonal, 9F13E4) HRP Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | WB, IHC, ELISA |
| Clonality | Monoclonal 9F13E4 |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Mouse |
| Uniprot ID | O75439 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human MPPB/PMPCB recombinant protein (Position: E23-Q479). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Mouse IgG2b |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | HRP |
| Suggested Dilutions | Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MPPB/PMPCB Antibody (monoclonal, 9F13E4) - HRP

For Research Use Only. Not for use in diagnostic procedures.