

Anti-PPWD1 Mouse Monoclonal Antibody [Clone ID: OTI1B2]

Catalog Number: M11798


Overview

Product Name	Anti-PPWD1 Mouse Monoclonal Antibody [Clone ID: OTI1B2]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio PPWD1 mouse monoclonal antibody, clone OTI1B2 (formerly 1B2). Catalog# M11798. Tested in IF, IHC, WB. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, WB
Clonality	Monoclonal OTI1B2
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q96BP3


Technical Details

Immunogen	Full length human recombinant protein of human PPWD1 (NP_056157) produced in HEK293T cell.
Isotype	IgG1
Concentration	0.6 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IF 1:100 IHC: 1:150


Anti-PPWD1 Mouse Monoclonal Antibody [Clone ID: OTI1B2] (M11798) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PPWD1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PPWD1.




Immunohistochemical staining of paraffin-embedded Human bladder tissue within the normal limits using anti-PPWD1 mouse monoclonal antibody. (M11798)



Immunohistochemical staining of paraffin-embedded Human Kidney tissue within the normal limits using anti-PPWD1 mouse monoclonal antibody. (M11798)



Immunohistochemical staining of paraffin-embedded Human liver tissue within the normal limits using anti-PPWD1 mouse monoclonal antibody. (M11798)



Anti-PPWD1 mouse monoclonal antibody (M11798) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY PPWD1.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PPWD1 Mouse Monoclonal Antibody [Clone ID: OTI1B2]

For Research Use Only. Not for use in diagnostic procedures.