

Anti-U2AF1L3 (U2AF1L4) Mouse Monoclonal Antibody [Clone ID: OTI8F10]

Catalog Number: M11846


Overview

Product Name	Anti-U2AF1L3 (U2AF1L4) Mouse Monoclonal Antibody [Clone ID: OTI8F10]
Reactive Species	Human
Description	Boster Bio U2AF1L4 mouse monoclonal antibody, clone OTI8F10. Catalog# M11846. Tested in WB. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal OTI8F10
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q8WU68

Technical Details

Immunogen	Full length human recombinant protein of human U2AF1L4 (NP_659424) produced in E.coli.
Isotype	IgG2a
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000

Anti-U2AF1L3 (U2AF1L4) Mouse Monoclonal Antibody [Clone ID: OTI8F10] (M11846) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY U2AF1L4 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-U2AF1L4 (1:2000).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-U2AF1L3 (U2AF1L4) Mouse Monoclonal Antibody [Clone ID: OTI8F10]

For Research Use Only. Not for use in diagnostic procedures.