

Anti-Meis homeobox 3 (MEIS3) Mouse Monoclonal Antibody [Clone ID: OTI2E12]

Catalog Number: M11850-1


Overview

Product Name	Anti-Meis homeobox 3 (MEIS3) Mouse Monoclonal Antibody [Clone ID: OTI2E12]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio MEIS3 (Meis homeobox 3) mouse monoclonal antibody, clone OTI2E12 (formerly 2E12). Catalog# M11850-1. Tested in FC, IF, IHC, WB. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IF, IHC, WB
Clonality	Monoclonal OTI2E12
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q99687


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 1-261 of human MEIS3 (NP_064545) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IHC 1:150 IF 1:100 Flow Cytometry 1:100


Anti-Meis homeobox 3 (MEIS3) Mouse Monoclonal Antibody [Clone ID: OTI2E12] (M11850-1) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MEIS3 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MEIS3.



Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-MEIS3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)





Immunohistochemical staining of paraffin-embedded Human pancreas tissue within the normal limits using anti-MEIS3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-MEIS3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)

HEK293T cells transfected with either MEIS3 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-MEIS3 antibody (M11850-1)



Anti-MEIS3 mouse monoclonal antibody (M11850-1) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY MEIS3.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Meis homeobox 3 (MEIS3) Mouse Monoclonal Antibody [Clone ID: OT12E12]

For Research Use Only. Not for use in diagnostic procedures.