

Anti-MPP3 Mouse Monoclonal Antibody [Clone ID: OTI8F1]

Catalog Number: M11924


Overview

Product Name	Anti-MPP3 Mouse Monoclonal Antibody [Clone ID: OTI8F1]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio MPP3 mouse monoclonal antibody, clone OTI8F1 (formerly 8F1). Catalog# M11924. Tested in FC, IF, WB. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IF, WB
Clonality	Monoclonal OTI8F1
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q13368


Technical Details

Immunogen	Full length human recombinant protein of human MPP3 (NP_001923) produced in HEK293T cell.
Isotype	IgG2b
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IF 1:100 Flow Cytometry 1:100


Anti-MPP3 Mouse Monoclonal Antibody [Clone ID: OTI8F1] (M11924) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MPP3 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MPP3.



Flow cytometric Analysis of Jurkat cells



Anti-MPP3 mouse monoclonal antibody (M11924) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY MPP3.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MPP3 Mouse Monoclonal Antibody [Clone ID: OTI8F1]

For Research Use Only. Not for use in diagnostic procedures.