

## Anti-ZCCHC17 Mouse Monoclonal Antibody [Clone ID: OTI1F11]

Catalog Number: M12008


### Overview

Product Name	Anti-ZCCHC17 Mouse Monoclonal Antibody [Clone ID: OTI1F11]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio ZCCHC17 mouse monoclonal antibody, clone OTI1F11. Catalog# M12008. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI1F11
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9NP64


### Technical Details

Immunogen	Full length human recombinant protein of human ZCCHC17 (NP_057589) produced in E.coli.
Isotype	IgG2a
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:

## Anti-ZCCHC17 Mouse Monoclonal Antibody [Clone ID: OTI1F11] (M12008) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZCCHC17 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZCCHC17 (1:2000).



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-ZCCHC17 monoclonal antibody (1:500).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZCCHC17 Mouse Monoclonal Antibody [Clone ID: OTI1F11]