

## Anti-GMFG Antibody(Ascites)

Catalog Number: M12403


### Overview

Product Name	Anti-GMFG Antibody(Ascites)
Reactive Species	Human
Description	Boster Bio Anti-GMFG Antibody(Ascites) (Catalog # M12403). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 789CT20.3.1.1
Formulation	Mouse monoclonal antibody supplied in crude ascites with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	O60234

### Technical Details

Immunogen	Purified His-tagged GMFG protein was used to produced this monoclonal antibody.
Isotype	IgG1
Suggested Dilutions	WB: 1:40000

## Anti-GMFG Antibody(Ascites) (M12403) Images



GMFG Antibody (N-term) western blot analysis in Jurkat cell line lysates (35ug/lane). This demonstrates the GMFG antibody detected the GMFG protein (arrow).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GMFG Antibody(Ascites)

For Research Use Only. Not for use in diagnostic procedures.