

Anti-FATE1 Mouse Monoclonal Antibody [Clone ID: OTI6F10]

Catalog Number: M12412-2


Overview

| | |
|----------------------|--|
| Product Name | Anti-FATE1 Mouse Monoclonal Antibody [Clone ID: OTI6F10] |
| Reactive Species | Dog, Human, Monkey, Mouse, Rat |
| Description | Boster Bio FATE1 mouse monoclonal antibody, clone OTI6F10 (formerly 6F10). Catalog# M12412-2. Tested in IF, WB. This antibody reacts with Human, Monkey, Mouse, Rat, Dog. |
| Application | IF, WB |
| Clonality | Monoclonal OTI6F10 |
| Formulation | PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C as received. |
| Host | Mouse |
| Uniprot ID | Q969F0 |


Technical Details

| | |
|---------------------|---|
| Immunogen | Full length human recombinant protein of human FATE1 (NP_001756) produced in HEK293T cell. |
| Isotype | IgG1 |
| Concentration | 1 mg/ml |
| Purification | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Suggested Dilutions | WB 1:200~4000 IF 1:100 |


Anti-FATE1 Mouse Monoclonal Antibody [Clone ID: OTI6F10] (M12412-2) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FATE1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FATE1.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-FATE1 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).



Anti-FATE1 mouse monoclonal antibody (M12412-2) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY FATE1.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FATE1 Mouse Monoclonal Antibody [Clone ID: OTI6F10]

For Research Use Only. Not for use in diagnostic procedures.