

## Anti-FATE1 Mouse Monoclonal Antibody [Clone ID: OTI6F10]

Catalog Number: M12412-2


### Overview

Product Name	Anti-FATE1 Mouse Monoclonal Antibody [Clone ID: OTI6F10]
Reactive Species	Dog, Human, Monkey, Mouse, Rat
Description	Boster Bio FATE1 mouse monoclonal antibody, clone OTI6F10 (formerly 6F10). Catalog# M12412-2. Tested in IF, WB. This antibody reacts with Human, Monkey, Mouse, Rat, Dog.
Application	IF, WB
Clonality	Monoclonal OTI6F10
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q969F0


### Technical Details

Immunogen	Full length human recombinant protein of human FATE1 (NP_001756) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:


## Anti-FATE1 Mouse Monoclonal Antibody [Clone ID: OTI6F10] (M12412-2) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FATE1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FATE1.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-FATE1 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).



Anti-FATE1 mouse monoclonal antibody (M12412-2) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY FATE1.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FATE1 Mouse Monoclonal Antibody [Clone ID: OTI6F10]