

Anti-Histone H3 (mono methyl R2) HIST1H3A Rabbit Monoclonal Antibody

Catalog Number: M12477-3


Overview

Product Name	Anti-Histone H3 (mono methyl R2) HIST1H3A Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Histone H3 (mono methyl R2) HIST1H3A Rabbit Monoclonal Antibody catalog # M12477-3. Tested in WB, ICC/IF applications. This antibody reacts with Human, Mouse.
Application	IF, ICC, WB
Clonality	Monoclonal DFB-8
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P68431


Technical Details

Immunogen	A synthesized peptide derived from human Histone H3 (mono methyl R2)
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 ICC/IF 1:50-200


Anti-Histone H3 (mono methyl R2) HIST1H3A Rabbit Monoclonal Antibody (M12477-3) Images



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:150 dilution.



Western blot analysis of Histone H3 (mono methyl R2) expression in HeLa cell lysate.

1 Publications Citing This Product

1. PubMed ID: 24503896, Role of endoplasmic reticulum stress in apoptosis of differentiated mouse podocytes induced by high glucose

Visit bosterbio.com/anti-histone-h3-mono-methyl-r2-rabbit-monoclonal-antibody-m12477-3-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Histone H3 (mono methyl R2) HIST1H3A Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.