

Anti-Histone H3 (mono+di+tri methyl K14) HIST1H3A Rabbit Monoclonal Antibody

Catalog Number: M12477-9


Overview

| | |
|----------------------|--|
| Product Name | Anti-Histone H3 (mono+di+tri methyl K14) HIST1H3A Rabbit Monoclonal Antibody |
| Reactive Species | Human, Mouse |
| Description | Boster Bio Anti-Histone H3 (mono+di+tri methyl K14) HIST1H3A Rabbit Monoclonal Antibody catalog # M12477-9. Tested in WB, ICC/IF applications. This antibody reacts with Human, Mouse. |
| Application | IF, ICC, WB |
| Clonality | Monoclonal I-8 |
| Formulation | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P68431/P84243/Q16695 |

Technical Details

| | |
|---------------------|--|
| Immunogen | A synthesized peptide derived from human Histone H3 (mono+di+tri methyl K14) |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5mg/ml |
| Purification | Affinity-chromatography |
| Suggested Dilutions | WB 1:500-2000 ICC/IF 1:50-200 |

Anti-Histone H3 (mono+di+tri methyl K14) HIST1H3A Rabbit Monoclonal Antibody (M12477-9) Images



Western blot analysis of Histone H3 (mono+di+tri methyl K14) expression in HeLa cell lysate.

1 Publications Citing This Product

1. PubMed ID: 24503896, Role of endoplasmic reticulum stress in apoptosis of differentiated mouse podocytes induced by high glucose

Visit bosterbio.com/anti-histone-h3-mono-di-tri-methyl-k14-rabbit-monoclonal-antibody-m12477-9-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Histone H3 (mono+di+tri methyl K14) HIST1H3A Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.