

Anti-COPS7B Antibody (Monoclonal, 34C04)

Catalog Number: M12717-1

About COPS7B

Involved in protein deneddylation. Located in nucleoplasm. Part of COP9 signalosome.

Overview

Product Name	Anti-COPS7B Antibody (Monoclonal, 34C04)
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-COPS7B Antibody (Monoclonal, 34C04) catalog # M12717-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal 34C04
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H9Q2

Technical Details

Immunogen	Synthetic peptide within human CSN7b aa 217-262.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-COPS7B Antibody (Monoclonal, 34C04)

For Research Use Only. Not for use in diagnostic procedures.