

## Anti-PAAF1 Mouse Monoclonal Antibody [Clone ID: OTI1F11]

Catalog Number: M12843


### Overview

Product Name	Anti-PAAF1 Mouse Monoclonal Antibody [Clone ID: OTI1F11]
Reactive Species	Human, Mouse
Description	Boster Bio PAAF1 mouse monoclonal antibody, clone OTI1F11 (formerly 1F11). Catalog# M12843. Tested in IF, WB. This antibody reacts with Human, Mouse.
Application	IF, WB
Clonality	Monoclonal OTI1F11
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9BRP4


### Technical Details

Immunogen	Full length human recombinant protein of human PAAF1 (NP_079431) produced in HEK293T cell.
Isotype	IgG2b
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:400~4000 IF 1:100


## Anti-PAAF1 Mouse Monoclonal Antibody [Clone ID: OTI1F11] (M12843) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PAAF1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PAAF1.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-PAAF1 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).



Anti-PAAF1 mouse monoclonal antibody (M12843) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY PAAF1.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PAAF1 Mouse Monoclonal Antibody [Clone ID: OTI1F11]

For Research Use Only. Not for use in diagnostic procedures.