

Anti-PPAP2B Antibody (N-term)

Catalog Number: M12857

About PLPP3

Catalyzes the conversion of phosphatidic acid (PA) to diacylglycerol (DG). In addition it hydrolyzes lysophosphatidic acid (LPA), ceramide-1-phosphate (C-1-P) and sphingosine-1-phosphate (S-1-P). The relative catalytic efficiency is LPA = PA > C-1-P > S-1-P. May be involved in cell adhesion and in cell-cell interactions.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-PPAP2B Antibody (N-term)
Reactive Species	Human, Rat
Description	Boster Bio Anti-PPAP2B Antibody (N-term) (Catalog # M12857). Tested in WB application(s). This antibody reacts with Human, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O14495

Technical Details

Immunogen	This PPAP2B antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 3~37 amino acids from the N-terminal region of human PPAP2B.
Predicted Reactive Species	Bovine
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.

Suggested Dilutions	WB: 1:2000
---------------------	------------

Anti-PPAP2B Antibody (N-term) (M12857) Images



All lanes : Anti-PPAP2B Antibody (N-term) at 1:2000 dilution
Lane 1: SW620 whole cell lysates
Lane 2: human placenta lysates
Lane 3: human kidney lysates
Lane 4: human lung lysates
Lane 5: rat liver lysates
Lysates/proteins at 20 µg per lane. Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/10000 dilution
Predicted band size : 35 kDa
Blocking/Dilution buffer: 5% NFDM/TBST.

Anti-PPAP2B Antibody (N-term)

For Research Use Only. Not for use in diagnostic procedures.