

Anti-PHYHD1 Rabbit Monoclonal Antibody

Catalog Number: M12964-1


Overview

| | |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Anti-PHYHD1 Rabbit Monoclonal Antibody |
| Reactive Species | Human |
| Description | This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of PHYHD1 protein. Validated for research with WB,IF/ICC,Flow Cyt applications. Anti PHYHD1 is reactive for Human samples |
| Application | Flow Cytometry, IF, ICC, WB |
| Clonality | Monoclonal |
| Formulation | PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C/1 year |
| Host | Rabbit |
| Uniprot ID | Q5SRE7 |

Technical Details

| | |
|---------------------|----------------------------------------------------------|
| Isotype | IgG, Kappa |
| Purification | Protein A |
| Suggested Dilutions | WB 1:1000-1:5000 ICC/IF 1:100-1:300 FC 1:100-1:300 |

Anti-PHYHD1 Rabbit Monoclonal Antibody (M12964-1) Images



Western Blot analysis of A549 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-PHYHD1 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PHYHD1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.