

Anti-GPC2 / Glypican 2 Reference Antibody (Nih Patent Anti-Glypican-2)

Catalog Number: M13075


Overview

Product Name	Anti-GPC2 / Glypican 2 Reference Antibody (Nih Patent Anti-Glypican-2)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-GPC2 / Glypican 2 Reference Antibody (Nih Patent Anti-Glypican-2) (Catalog # M13075). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.992EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q8N158


Technical Details

Isotype	IgG1
Form	Liquid

Anti-GPC2 / Glypican 2 Reference Antibody (Nih Patent Anti-Glypican-2) (M13075) Images



Anti-GPC2/Glypican 2 Reference Antibody (Nih Patent Anti-Glypican-2) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-GPC2/Glypican 2 Reference Antibody (Nih Patent Anti-Glypican-2) is more than 95%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GPC2 / Glypican 2 Reference Antibody (Nih Patent Anti-Glypican-2)

For Research Use Only. Not for use in diagnostic procedures.