

Anti-UBE2D4 Mouse Monoclonal Antibody [Clone ID: OTI1B1]

Catalog Number: M13194


Overview

Product Name	Anti-UBE2D4 Mouse Monoclonal Antibody [Clone ID: OTI1B1]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio UBE2D4 mouse monoclonal antibody, clone OTI1B1. Catalog# M13194. Tested in IHC, WB. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal OTI1B1
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9Y2X8


Technical Details

Immunogen	Full length human recombinant protein of human UBE2D4 (NP_057067) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500~2000 IHC 1:150


Anti-UBE2D4 Mouse Monoclonal Antibody [Clone ID: OTI1B1] (M13194) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY UBE2D4 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UBE2D4 (1:2000).



Immunohistochemical staining of paraffin-embedded Human adult heart tissue within the normal limits using anti-UBE2D4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)



Immunohistochemical staining of paraffin-embedded Human embryonic cerebellum within the normal limits using anti-UBE2D4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UBE2D4 Mouse Monoclonal Antibody [Clone ID: OTI1B1]

For Research Use Only. Not for use in diagnostic procedures.