

Anti-PCTAIRE2 (CDK17) Mouse Monoclonal Antibody [Clone ID: OTI2E11]

Catalog Number: M13416


Overview

Product Name	Anti-PCTAIRE2 (CDK17) Mouse Monoclonal Antibody [Clone ID: OTI2E11]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio CDK17 mouse monoclonal antibody, clone OTI2E11. Catalog# M13416. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI2E11
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q00537

Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 6-132 of human CDK17 (NP_001163935) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)

Anti-PCTAIRE2 (CDK17) Mouse Monoclonal Antibody [Clone ID: OTI2E11] (M13416) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CDK17 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CDK17.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PCTAIRE2 (CDK17) Mouse Monoclonal Antibody [Clone ID: OTI2E11]