

Anti-Golgin 97/GOLGA1 Antibody Picoband® (monoclonal, 8E4H1) HRP Conjugated

Catalog Number: M13524-1-HRP

About GOLGA1

Golgin subfamily A member 1 is a protein that in humans is encoded by the GOLGA1 gene. The Golgi apparatus, which participates in glycosylation and transport of proteins and lipids in the secretory pathway, consists of a series of stacked cisternae (flattened membrane sacs). Interactions between the Golgi and microtubules are thought to be important for the reorganization of the Golgi after it fragments during mitosis. This gene encodes one of the golgins, a family of proteins localized to the Golgi. This encoded protein is associated with Sjogren's syndrome.

Overview

Product Name	Anti-Golgin 97/GOLGA1 Antibody Picoband® (monoclonal, 8E4H1) HRP Conjugated
Reactive Species	Human
Clonality	Monoclonal 8E4H1
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Mouse
Uniprot ID	Q92805

Technical Details

Immunogen	E.coli-derived human Golgin 97/GOLGA1 recombinant protein (Position: M1-K752).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	IgG1
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Golgin 97/GOLGA1 Antibody (monoclonal, 8E4H1) - HRP

For Research Use Only. Not for use in diagnostic procedures.