

## Anti-PANK3 Mouse Monoclonal Antibody [Clone ID: OTI3C4]

Catalog Number: M13559


### Overview

Product Name	Anti-PANK3 Mouse Monoclonal Antibody [Clone ID: OTI3C4]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio PANK3 mouse monoclonal antibody, clone OTI3C4 (formerly 3C4). Catalog# M13559. Tested in FC, WB. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, WB
Clonality	Monoclonal OTI3C4
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9H999


### Technical Details

Immunogen	Full length human recombinant protein of human PANK3 (NP_078870) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500 Flow Cytometry 1:100


## Anti-PANK3 Mouse Monoclonal Antibody [Clone ID: OTI3C4] (M13559) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PANK3 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PANK3.



HEK293T cells transfected with either PANK3 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-PANK3 antibody (M13559)



Flow cytometric Analysis of HeLa cells



Flow cytometric Analysis of Jurkat cells

**Submit a product review to [Biocompare.com](https://www.biocompare.com)**

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PANK3 Mouse Monoclonal Antibody [Clone ID: OTI3C4]

For Research Use Only. Not for use in diagnostic procedures.