

Anti-DNAAF10 Mouse Monoclonal Antibody [Clone ID: OTI4B5]

Catalog Number: M13880


Overview

Product Name	Anti-DNAAF10 Mouse Monoclonal Antibody [Clone ID: OTI4B5]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio WDR92 mouse monoclonal antibody, clone OTI4B5. Catalog# M13880. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI4B5
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q96MX6


Technical Details

Immunogen	Full length human recombinant protein of human WDR92 (NP_612467) produced in HEK293T cell.
Isotype	IgG2b
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500~2000

Anti-DNAAF10 Mouse Monoclonal Antibody [Clone ID: OTI4B5] (M13880) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY WDR92 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-WDR92 (1:2000).



Western blot analysis of extracts (35ug) from 2 different cell lines and human kidney tissue lysate by using anti-WDR92 monoclonal antibody (1:500).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DNAAF10 Mouse Monoclonal Antibody [Clone ID: OTI4B5]

For Research Use Only. Not for use in diagnostic procedures.