

Anti-Heat Shock Factor 2 Binding Protein (HSF2BP) Mouse Monoclonal Antibody [Clone ID: OTI2A1]

Catalog Number: M14085


Overview

| | |
|----------------------|--|
| Product Name | Anti-Heat Shock Factor 2 Binding Protein (HSF2BP) Mouse Monoclonal Antibody [Clone ID: OTI2A1] |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio HSF2BP mouse monoclonal antibody, clone OTI2A1. Catalog# M14085. Tested in WB. This antibody reacts with Human, Mouse, Rat. |
| Application | WB |
| Clonality | Monoclonal OTI2A1 |
| Formulation | PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C as received. |
| Host | Mouse |
| Uniprot ID | O75031 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Human recombinant protein fragment corresponding to amino acids 1-181 of human HSF2BP (NP_008962) produced in E.coli. |
| Isotype | IgG1 |
| Concentration | 1 mg/ml |
| Purification | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Suggested Dilutions | WB 1:500 |

Anti-Heat Shock Factor 2 Binding Protein (HSF2BP) Mouse Monoclonal Antibody [Clone ID: OTI2A1] (M14085) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HSF2BP (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HSF2BP (1:500).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Heat Shock Factor 2 Binding Protein (HSF2BP) Mouse Monoclonal Antibody [Clone ID: OTI2A1]

For Research Use Only. Not for use in diagnostic procedures.