

Anti-KBTBD7 Mouse Monoclonal Antibody [Clone ID: OTI1G3]

Catalog Number: M14135


Overview

Product Name	Anti-KBTBD7 Mouse Monoclonal Antibody [Clone ID: OTI1G3]
Reactive Species	Dog, Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-KBTBD7 mouse monoclonal antibody, clone OTI1G3 (formerly 1G3). Catalog# M14135. Tested in IF, WB. This antibody reacts with Human, Monkey, Mouse, Rat, Dog.
Application	IF, WB
Clonality	Monoclonal OTI1G3
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q8WVZ9


Technical Details

Immunogen	Full length human recombinant protein of human KBTBD7 (NP_115514) produced in HEK293T cell.
Isotype	IgG2b
Concentration	0.84 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IF 1:100


Anti-KBTBD7 Mouse Monoclonal Antibody [Clone ID: OTI1G3] (M14135) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY KBTBD7 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-KBTBD7.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-KBTBD7 monoclonal antibody.



Anti-KBTBD7 mouse monoclonal antibody (M14135) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY KBTBD7.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KBTBD7 Mouse Monoclonal Antibody [Clone ID: OTI1G3]

For Research Use Only. Not for use in diagnostic procedures.