

## Anti-PLEK2 Mouse Monoclonal Antibody [Clone ID: OTI5D3]

Catalog Number: M14268


### Overview

Product Name	Anti-PLEK2 Mouse Monoclonal Antibody [Clone ID: OTI5D3]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio PLEK2 mouse monoclonal antibody, clone OTI5D3 (formerly 5D3). Catalog# M14268. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI5D3
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9NYT0


### Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 132-353 of human PLEK2 (NP_057529) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000

## Anti-PLEK2 Mouse Monoclonal Antibody [Clone ID: OTI5D3] (M14268) Images



Western blot analysis of HT29 cell lysate (35ug) by using anti-PLEK2 monoclonal antibody. Dilution: 1:500



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PLEK2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PLEK2.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PLEK2 Mouse Monoclonal Antibody [Clone ID: OTI5D3]

For Research Use Only. Not for use in diagnostic procedures.