

Anti-TPSG1 Mouse Monoclonal Antibody [Clone ID: OTI2A10]

Catalog Number: M14320-1


Overview

Product Name	Anti-TPSG1 Mouse Monoclonal Antibody [Clone ID: OTI2A10]
Reactive Species	Human
Description	Boster Bio TPSG1 mouse monoclonal antibody, clone OTI2A10 (formerly 2A10). Catalog# M14320-1. Tested in FC, IF, WB. This antibody reacts with Human.
Application	Flow Cytometry, IF, WB
Clonality	Monoclonal OTI2A10
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9NRR2


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 20-283 of human TPSG1 (NP_036599) produced in E.coli.
Isotype	IgG1
Concentration	0.52 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IF 1:100 Flow Cytometry 1:100


Anti-TPSG1 Mouse Monoclonal Antibody [Clone ID: OT12A10] (M14320-1) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TPSG1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TPSG1.



HEK293T cells transfected with either TPSG1 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-TPSG1 antibody (M14320-1)



Anti-TPSG1 mouse monoclonal antibody (M14320-1) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY TPSG1.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TPSG1 Mouse Monoclonal Antibody [Clone ID: OT12A10]

For Research Use Only. Not for use in diagnostic procedures.