

## Anti-TOMM40L Mouse Monoclonal Antibody [Clone ID: OTI5B1]

Catalog Number: M14446


### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-TOMM40L Mouse Monoclonal Antibody [Clone ID: OTI5B1]   |
| Reactive Species     | Human, Mouse, Rat   |
| Description          | Boster Bio TOMM40L mouse monoclonal antibody, clone OTI5B1. Catalog# M14446. Tested in WB. This antibody reacts with Human, Mouse, Rat. |
| Application          | WB  |
| Clonality            | Monoclonal OTI5B1   |
| Formulation          | PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.  |
| Storage Instructions | Store at -20°C as received.   |
| Host                 | Mouse   |
| Uniprot ID           | Q969M1  |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Full length human recombinant protein of human TOMM40L (NP_115550) produced in E.coli.  |
| Isotype             | IgG2b   |
| Concentration       | 1 mg/ml   |
| Purification        | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)   |
| Suggested Dilutions | Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.<br>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.<br>Some PubMed article(s) citing the expression level of this target are as follows:<br><br>For protocols, please visit <a href="https://www.bosterbio.com/protocol-and-troubleshooting/">https://www.bosterbio.com/protocol-and-troubleshooting/</a> |

## Anti-TOMM40L Mouse Monoclonal Antibody [Clone ID: OTI5B1] (M14446) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TOMM40L (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TOMM40L (1:500).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TOMM40L Mouse Monoclonal Antibody [Clone ID: OTI5B1]