

Anti-Kelch-like protein 25 KLHL25 Monoclonal Antibody

Catalog Number: M14745


Overview

Product Name	Anti-Kelch-like protein 25 KLHL25 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Kelch-like protein 25 KLHL25 Monoclonal Antibody catalog # M14745. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 1B7
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9H0H3

Technical Details

Immunogen	Purified recombinant fragment of human KLHL25 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-Kelch-like protein 25 KLHL25 Monoclonal Antibody (M14745) Images



Western Blot analysis using KLHL25 Monoclonal Antibody against KLHL25-hlgGfc transfected HEK293 cell.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Kelch-like protein 25 KLHL25 Monoclonal Antibody
For Research Use Only. Not for use in diagnostic procedures.