

Anti-ZNF213 Mouse Monoclonal Antibody [Clone ID: OTI1A12]

Catalog Number: M15093


Overview

Product Name	Anti-ZNF213 Mouse Monoclonal Antibody [Clone ID: OTI1A12]
Reactive Species	Human
Description	Boster Bio ZNF213 mouse monoclonal antibody, clone OTI1A12. Catalog# M15093. Tested in WB. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal OTI1A12
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	O14771

Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 1-292 of human ZNF213 (NP_004211) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:200

Anti-ZNF213 Mouse Monoclonal Antibody [Clone ID: OTI1A12] (M15093) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZNF213 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF213 (1:200).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZNF213 Mouse Monoclonal Antibody [Clone ID: OTI1A12]

For Research Use Only. Not for use in diagnostic procedures.