

Anti-CLEC1 (CLEC1A) Mouse Monoclonal Antibody [Clone ID: OTI1H8]

Catalog Number: M15410


Overview

Product Name	Anti-CLEC1 (CLEC1A) Mouse Monoclonal Antibody [Clone ID: OTI1H8]
Reactive Species	Human
Description	Boster Bio CLEC1A mouse monoclonal antibody, clone OTI1H8. Catalog# M15410. Tested in WB. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal OTI1H8
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q8NC01

Technical Details

Immunogen	Full length human recombinant protein of human CLEC1A (NP_057595) produced in HEK293T cell.
Isotype	IgG2b
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000

Anti-CLEC1 (CLEC1A) Mouse Monoclonal Antibody [Clone ID: OTI1H8] (M15410) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CLEC1A (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CLEC1A (1:2000).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CLEC1 (CLEC1A) Mouse Monoclonal Antibody [Clone ID: OTI1H8]

For Research Use Only. Not for use in diagnostic procedures.