

Anti-C20orf30 (TMEM230) Mouse Monoclonal Antibody [Clone ID: OTI2C8]

Catalog Number: M15980


Overview

Product Name	Anti-C20orf30 (TMEM230) Mouse Monoclonal Antibody [Clone ID: OTI2C8]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio C20orf30 mouse monoclonal antibody, clone OTI2C8 (formerly 2C8). Catalog# M15980. Tested in FC, IF, WB. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IF, WB
Clonality	Monoclonal OTI2C8
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q96A57


Technical Details

Immunogen	Full length human recombinant protein of human C20orf30 (NP_054864) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:


Anti-C20orf30 (TMEM230) Mouse Monoclonal Antibody [Clone ID: OTI2C8] (M15980) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY C20orf30 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-C20orf30.



HEK293T cells transfected with either TMEM230 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-C20orf30 antibody (M15980)



Anti-C20orf30 mouse monoclonal antibody (M15980) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY C20orf30.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C20orf30 (TMEM230) Mouse Monoclonal Antibody [Clone ID: OTI2C8]