

Anti-ZNF114 Mouse Monoclonal Antibody [Clone ID: OTI6F4]

Catalog Number: M16543


Overview

| | |
|----------------------|--|
| Product Name | Anti-ZNF114 Mouse Monoclonal Antibody [Clone ID: OTI6F4] |
| Reactive Species | Human |
| Description | Boster Bio ZNF114 mouse monoclonal antibody, clone OTI6F4 (formerly 6F4). Catalog# M16543. Tested in WB. This antibody reacts with Human. |
| Application | WB |
| Clonality | Monoclonal OTI6F4 |
| Formulation | PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C as received. |
| Host | Mouse |
| Uniprot ID | Q8NC26 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Full length human recombinant protein of human ZNF114 (NP_705836) produced in E.coli. |
| Isotype | IgG1 |
| Concentration | 1 mg/ml |
| Purification | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Suggested Dilutions | WB 1:2000 |

Anti-ZNF114 Mouse Monoclonal Antibody [Clone ID: OTI6F4] (M16543) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZNF114 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF114 (1:2000).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZNF114 Mouse Monoclonal Antibody [Clone ID: OTI6F4]

For Research Use Only. Not for use in diagnostic procedures.