

Anti-UBXD3 (UBXN10) Mouse Monoclonal Antibody [Clone ID: OTI1H5]

Catalog Number: M16719


Overview

Product Name	Anti-UBXD3 (UBXN10) Mouse Monoclonal Antibody [Clone ID: OTI1H5]
Reactive Species	Human
Description	Boster Bio Anti-UBXN10 mouse monoclonal antibody, clone OTI1H5 (formerly 1H5). Catalog# M16719. Tested in FC, IF, WB. This antibody reacts with Human.
Application	Flow Cytometry, IF, WB
Clonality	Monoclonal OTI1H5
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q96LJ8


Technical Details

Immunogen	Full length human recombinant protein of human UBXN10 (NP_689589) produced in HEK293T cell.
Isotype	IgG1
Concentration	0.45 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500 IF 1:100 Flow Cytometry 1:100


Anti-UBXD3 (UBXN10) Mouse Monoclonal Antibody [Clone ID: OT11H5] (M16719) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY UBXN10 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UBXN10.



HEK293T cells transfected with either UBXN10 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-UBXN10 antibody (M16719)



Anti-UBXN10 mouse monoclonal antibody (M16719) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY UBXN10.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UBXD3 (UBXN10) Mouse Monoclonal Antibody [Clone ID: OT11H5]

For Research Use Only. Not for use in diagnostic procedures.