

## Anti-ZMAT4 Mouse Monoclonal Antibody [Clone ID: OTI9E11]

Catalog Number: M17345


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-ZMAT4 Mouse Monoclonal Antibody [Clone ID: OTI9E11]   |
| Reactive Species     | Human, Mouse, Rat  |
| Description          | Boster Bio ZMAT4 mouse monoclonal antibody, clone OTI9E11. Catalog# M17345. Tested in WB. This antibody reacts with Human, Mouse, Rat.   |
| Application          | WB   |
| Clonality            | Monoclonal OTI9E11   |
| Formulation          | PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C as received.  |
| Host                 | Mouse  |
| Uniprot ID           | Q9H898   |


### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Full length human recombinant protein of human ZMAT4 (NP_078921) produced in E.coli.                      |
| Isotype             | IgG1  |
| Concentration       | 1 mg/ml   |
| Purification        | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Suggested Dilutions | WB 1:500~2000   |

## Anti-ZMAT4 Mouse Monoclonal Antibody [Clone ID: OTI9E11] (M17345) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZMAT4 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZMAT4 (1:2000).



Western blot analysis of extracts (35ug) from HELA cell line by using anti-ZMAT4 monoclonal antibody (1:500).

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZMAT4 Mouse Monoclonal Antibody [Clone ID: OTI9E11]

For Research Use Only. Not for use in diagnostic procedures.