

## Anti-AVPI1 Mouse Monoclonal Antibody [Clone ID: OTI1G1]

Catalog Number: M18131


### Overview

Product Name	Anti-AVPI1 Mouse Monoclonal Antibody [Clone ID: OTI1G1]
Reactive Species	Human
Description	Boster Bio AVPI1 mouse monoclonal antibody, clone OTI1G1 (formerly 1G1). Catalog# M18131. Tested in IHC, WB. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal OTI1G1
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q5T686


### Technical Details

Immunogen	Full length human recombinant protein of human AVPI1 (NP_068378) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:

## Anti-AVPI1 Mouse Monoclonal Antibody [Clone ID: OTI1G1] (M18131) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY AVPI1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AVPI1.



Immunohistochemical staining of paraffin-embedded Human colon tissue within the normal limits using anti-AVPI1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AVPI1 Mouse Monoclonal Antibody [Clone ID: OTI1G1]