

## Anti-TBC1D28 Mouse Monoclonal Antibody [Clone ID: OTI4A9]

Catalog Number: M19702


### Overview

Product Name	Anti-TBC1D28 Mouse Monoclonal Antibody [Clone ID: OTI4A9]
Reactive Species	Human
Description	Boster Bio TBC1D28 mouse monoclonal antibody, clone OTI4A9 (formerly 4A9). Catalog# M19702. Tested in IHC, WB. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal OTI4A9
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q2M2D7


### Technical Details

Immunogen	Full length human recombinant protein of human TBC1D28 (NP_001034486) produced in E.coli.
Isotype	IgG2a
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IHC 1:150


## Anti-TBC1D28 Mouse Monoclonal Antibody [Clone ID: OTI4A9] (M19702) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TBC1D28 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TBC1D28 (1:2000).



Immunohistochemical staining of paraffin-embedded Human adult heart tissue using anti-TBC1D28 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)



Immunohistochemical staining of paraffin-embedded Human adult brain tissue using anti-TBC1D28 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TBC1D28 Mouse Monoclonal Antibody [Clone ID: OTI4A9]

For Research Use Only. Not for use in diagnostic procedures.