

Anti-GRK3 Monoclonal Antibody

Catalog Number: M20057


Overview

Product Name	Anti-GRK3 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-GRK3 Monoclonal Antibody catalog # M20057. Tested in WB, IHC, ICC/IF applications. This antibody reacts with Human.
Application	IF, IHC, ICC, WB
Clonality	Monoclonal AEBE-7
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P35626

Technical Details

Immunogen	A synthesized peptide derived from human GRK3 Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200

Anti-GRK3 Monoclonal Antibody (M20057) Images



Western blot analysis of GRK3 expression in Ramos cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GRK3 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.