

## Anti-hCD45R Antibody (Monoclonal, BRA-11G)

Catalog Number: M30397

### About PTPRC

By immunohistochemistry, it specifically recognizes a protein in melanocytes and melanomas. This monoclonal antibody reacts with junctional and blue nevus cells and variably with fetal and neonatal melanocytes. Intradermal nevi, normal adult melanocytes, and non-melanocytic cells are negative. It does not stain tumor cells of epithelial, lymphoid, glial, or mesenchymal origin. This monoclonal antibody labels formalin-fixed, paraffin-embedded melanomas and other tumors showing melanocytic differentiation.

### Overview


Product Name	Anti-hCD45R Antibody (Monoclonal, BRA-11G)
Reactive Species	Human
Description	Boster Bio Anti-hCD45R Antibody (Monoclonal, BRA-11G) (Catalog# M30397). Tested in FACS application(s). This antibody reacts with Human.
Conjugate	Biotin
Application	FACS
Clonality	Monoclonal BRA-11G
Formulation	100 µg in 200 µl PBS containing 0.05% stabilizing protein and 0.05% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store the antibody at 4°C, stable for 6 months. For long-term storage, store at -20°C. Avoid repeated freeze and thaw cycles.
Host	Mouse
Uniprot ID	P08575

### Technical Details

Immunogen	Non-T, non-B CALLA positive ALL cell line REH (Leucocyte Workshop IV and V)
Predicted Reactive Species	Pig, Rabbit
Cross Reactivity	Does not cross-react with primate, avian or amphibian GR.
Isotype	IgG1 Kappa
Form	Liquid

Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein G Chromatography
Suggested Dilutions	FACS: 0.5-1 $\mu\text{g}/10^6$ cells

## Anti-hCD45R Antibody (Monoclonal, BRA-11G) (M30397) Images



Cell Surface flow analysis of hCD45R in PBMC (Lymphocyte gated) using 0.5  $\mu$ g/ml

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-hCD45R Antibody (Monoclonal, BRA-11G)

For Research Use Only. Not for use in diagnostic procedures.